Beskrivelse

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
Description


Falling bodies apparatus or Karise cannon - It can be used to show that a throwing movement can be perceived as composed of a free fall and a horizontal movement. The movement is initiated by the release arm and the throwing mechanism fires one bullet and brings the other into free fall.

The two balls can be stored inside the base outside of use. Rotate the black plate to open the storage tray.

Technical data

- Dimensions: 180 x 170 x 50 mm.
- Ball diameter: 19 mm.
- Horizontal speed: Approx. 3.5 m / s.

Demostration use

The Karise cannon fires two bullets at the same time. Both spheres are affected vertically by gravity. One is simply released and falls vertically down. The second bullet is fired at a horizontal initial velocity; it falls

along a parabolic curve. When the balls hit the floor, it sounds like one almost "click".

Experiments

The "throwing motion" consists of both a free fall and a horizontal movement:

Demonstrate how the throwing motion can

be perceived as composed of free fall and a horizontal motion:

1. Place the Karise cannon on a horizontal table top so that the firing arm is outside the table edge as shown.
2. Pull back the firing lever until it is caught by the trigger.
3. Place a steel ball in both of the two small holes at the ends of the arm.
4. Press the shutter button. One bullet falls straight down; the other is moved horizontally.
5. Then, listen for the click it gives when the balls hit the floor. If both the Karise cannon and the floor are horizontal, the balls should hit the floor at simultaneously.

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© Frederiksen Scientific A/S, Ølgod.