Positiv feltretning

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The Tesla meter is used for measuring magnetic fields up to 2 T. The device consists of a loose probe and a measuring box with a digital reading of the result.


The device is powered by a 12 V adapter. A Ø 10 mm clamping pin is included, which can be mounted on the probe.

The magnetic field detects with a Hall element which is sensitive to magnetic fields perpendicular to the flat tip of the probe.

The probe is equipped with an 8-pin DIN connector, which also fits Pasco's data acquisition equipment (analog adapter PS-2158 or an interface UI-5000 and UI-5001 (550 or 850)).

Technical data

- Measuring range 100 mT – 1999 mT (Resolution 1 mT)
- Measuring range 0.0 – 199.9 mT (Resolution 0.1 mT)
- Automatic change of measurement range:
 - Upward is changed from 199.9 mT to 200 mT
 - Downward is changed from 100 mT to 99.9 mT
- Accuracy: 5%

Positive field direction

Calibration when used with data acquisition equipment: 5000 mV ~ 2000 mT.

Dimensions

- Magnetic field probe (tip): 80 x 8 x 2 mm
- Magnetic field probe (total): 215 x 35 x 25 mm
- Instrument: 158 x 108 x 56 mm.

Spare parts

355050 Mains adapter 12 V, with DC plug
406055 Magnetic field probe.

Experiment with the tesla meter

The following descriptions are intended as inspiration for experimental studies of magnetic fields and do not constitute complete practice instructions.

Magnetic field from a coil – varying current strength

Place a coil, e.g., 462520 (400 turns) next to the probe tip. Vary the current in 5 – 10 steps up to the maximum allowed for the coil (in the mentioned case 1 A).

Data processing

Plot data in a coordinate system. Describe how the magnetic field depends on the current.

The effect of an iron core.

Repeat the above with an iron core in the coil.

Magnetic field from a coil – varying number of turns.

Position the probe and a coil so that it is possible to replace the coil without changing their relative position.

Choose three – four coils (e.g., 462510, 462520, 462525 and 462530). Find the maximum current that all the coils can withstand and use this for the measurements. Measure the magnetic field for each coil.

Data processing

Plot the magnetic field as a function of the number of turns in a coordinate system. Describe the context.

The magnetic field near the end of a bar magnet

Place a bar magnet (e.g., 330510 - supplied as a pair, use only one) on a piece of graph paper so that it is easy to determine the distance between the end of the magnet and the Teslameter probe.

Measure the magnetic field along the extension of the magnet at a distance of 10 mm. Repeat in steps of 5 mm out to a distance corresponding to half the length of the magnet. Draw a graph of the measurement results.

Data processing (requires spreadsheets)

If you imagine that the bar magnet consists of two point-shaped magnetic poles, we have measured the field from one pole, while the other is "very far away". It would be natural to guess that the field decreases proportionally to the square of the distance.

- Calculate for all distances r the magnitude $1/r^2$ and plot the magnetic field B as a function of $1/r^2$.
- Are the points on a straight line?
- Try adding a constant distance r_0 to the distances and plot B as a function of $1/(r+r_0)^2$.
- Vary the size of r_0 until the points lie as close to a straight line as possible.
- Give an interpretation of the quantity r_0 .

Magnetic field from a dipole

Use a neodymium magnet (e.g., 331820 – supplied as a pair, use only one).

Place the magnet (possibly with a bit of tape) on a piece of graph paper so that it is easy to determine the distance between the end of the magnet and the Teslameter probe.

It must be measured along the extension of the center axis of the magnet.

First measure the magnetic field at a distance of 5 mm. Repeat in steps of 5 mm out to a distance of 70 mm.

Draw a graph of the measurement results.

Data processing (Requires spreadsheet)

As in the experiment above, an r_0 correction is introduced. Make a separate column with $r+r_0$.

Plot the measurement results in three double-logarithmic coordinate systems: One with all measurements, one with data from 5 to 20 mm, one with data from 30 to 70 mm.

In all three cases fit with a power function. Is it reasonable to describe the distance dependence as proportional to $1/r^2$ (with or without r_0)?