



The Biotechnology Education Company®

PO Box 341232
Bethesda, MD 20827

ph: 1.800.EDVOTEK (800.338.6835)
fax: 301.340.0582

email: info@edvotek.com
web: www.edvotek.com

IMPORTANT NOTICE – 5010 TetraSource™ Power Source

EDVOTEK® has designed the 5010 TetraSource™ Power Source to be the safest, most user-friendly, intelligent unit in its class. It features a pre-programmed safety threshold set at 60 Volts, which may require adjustment depending on your experimental conditions. If your experiment utilizes a low-draw electrophoresis buffer and an error message is displayed as you start your electrophoresis run, please follow these simple instructions below to adjust the safety volt threshold:

- While pressing the "**Δ VALUE**" button (upper right), press the "**ON**" button to turn on the power supply.
- Release the "**Δ VALUE**" button.
- Using the "**MODE**" button (upper left), scroll down to "**VOLT THRESHOLD**".
- Once at the "**VOLT THRESHOLD**", press the "**Δ VALUE**" to increase the voltage.

**Note – in most cases, the volt threshold will need to be increased. First, try increasing it to 120 volts.*

- Once the desired voltage is selected, use the "**MODE**" button to scroll down to "**EXIT**".
- Press the "**Δ VALUE**" button to exit the Setup Menu.

The 5010 TetraSource™ Power Supply should now be ready for your electrophoresis experiment. If additional adjustments are required, repeat the above steps to select a different volt threshold.

If you have any questions, please don't hesitate to contact us at 1.800.338.6835 or info@edvotek.com and we'll be happy to assist you.